KONICA MINOLTA

SONIMAGE HS1

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SONIMAGE HS1

**DICOM3.0
Conformance
Statement**

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内容	値
Implementation Class UID	1.2.392.200036.9107.700.1
Implementation Version Name	MG_DICOM V1.00

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3.4.2.3 DICOMDIR 属性

タグ	属性名	注釈
File-set Identification		
(0004,1130)	File-set ID	
Directory Information		
(0004,1200)	Offset of the First Directory Record of the Root Directory Entity	
(0004,1202)	Offset of the Last Directory Record of the Root Directory Entity	
(0004,1212)	File-set Consistency Flag	
(0004,1220)	Directory Record Sequence	
>(0004,1400)	Offset of the Next Directory Record	
>(0004,1410)	Record In-use Flag	
>(0004,1420)	Offset of Referenced Lower-Level Directory Entity	
>(0004,1430)	Directory Record Type	以下の列挙値を設定します。 PATIENT STUDY SERIES IMAGE
>(0004,1500)	Referenced File ID	Directory Record Type が IMAGE の時のみ出力します。
>(0004,1510)	Referenced SOP Class UID in File	Directory Record Type が IMAGE の時のみ出力します。
>(0004,1511)	Referenced SOP Instance UID in File	Directory Record Type が IMAGE の時のみ出力します。
>(0004,1512)	Referenced Transfer Syntax UID in File	Directory Record Type が IMAGE の時のみ出力します。
Patient Keys(Directory Record Type PATIENT)		
(0008,0005)	Specific Character Set	
(0010,0010)	Patient's Name	
(0010,0020)	Patient ID	
(0010,0030)	Patient's Birth Date	
(0010,0040)	Patient's Sex	
Study Keys(Directory Record Type STUDY)		
(0008,0005)	Specific Character Set	
(0008,0020)	Study Date	
(0008,0030)	Study Time	
(0008,0050)	Accession Number	
(0008,1030)	Study Description	空の値を設定します。
(0020,000D)	Study Instance UID	
(0020,0010)	Study ID	
Series Keys(Directory Record Type SERIES)		
(0008,0005)	Specific Character Set	
(0008,0060)	Modality	
(0020,000E)	Series Instance UID	
(0020,0011)	Series Number	
Image Keys(Directory Record Type IMAGE)		
(0008,0005)	Specific Character Set	
(0020,0013)	Instance Number	

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7.1.2.3 Vascular Patient Characteristics

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Age	DCM	121033	Subject Age	
Gender	DCM	121032	Subject Sex	

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Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
PSV	LN	11726-7	Peak Systolic Velocity	
Ved	LN	11653-3	End Diastolic Velocity	
S/D	LN	12144-2	Systolic to Diastolic Velocity Ratio	
RI	LN	12023-8	Resistivity Index	
PI	LN	12008-9	Pulsatility Index	
Vm_peak	LN	11692-1	Time averaged peak velocity	

7.1.2.5 Artery of Neck (Left Extracranial Arteries, Carotid Ratios)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Artery of Neck (Left Extracranial Arteries,Carotid Ratios)	SRT	T-45005	Artery of neck	Laterality: Left(SRT, G-A101)
Lt.CCA	SRT	T-45100	Common Carotid Artery	
Prox	7.1.2.4 を参照			Topographical Modifier: Proximal(SRT, G-A118) VesselBranch: Left(SRT, G-A101)
Prox ICA/CCA(S)	LN	33868-1	ICA/CCA velocity ratio	Topographical Modifier: Proximal(SRT, G-A118)
Prox ICA/CCA(D)	99KMHC	KHU-103-00001	ICA/CCA velocity ratio(Ved)	Topographical Modifier: Proximal(SRT, G-A118)
Mid	7.1.2.4 を参照			Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Left(SRT, G-A101)
Mid ICA/CCA(S)	LN	33868-1	ICA/CCA velocity ratio	Topographical Modifier: Mid-longitudinal(SRT, G-A188)
Mid ICA/CCA(D)	99KMHC	KHU-103-00001	ICA/CCA velocity ratio(Ved)	Topographical Modifier: Mid-longitudinal(SRT, G-A188)
Dist	7.1.2.4 を参照			Topographical Modifier: Distal(SRT, G-A119) VesselBranch: Left(SRT, G-A101)
Dist ICA/CCA(S)	LN	33868-1	ICA/CCA velocity ratio	Topographical Modifier: Distal(SRT, G-A119)
Dist ICA/CCA(D)	99KMHC	KHU-103-00001	ICA/CCA velocity ratio(Ved)	Topographical Modifier: Distal(SRT, G-A119)
Lt.BIF	SRT	T-45160	Carotid Bifurcation	
	7.1.2.4 を参照			VesselBranch: Left(SRT, G-A101)

Lt.Bulb	SRT	T-45170	Carotid Bulb	
	7.1.2.4 を参照			VesselBranch: Left(SRT, G-A101)
Lt.ECA	SRT	T-45200	External Carotid Artery	
	7.1.2.4 を参照			VesselBranch: Left(SRT, G-A101)
Lt.ICA	SRT	T-45300	Internal Carotid Artery	
Prox	7.1.2.4 を参照			Topographical Modifier: Proximal(SRT, G-A118) VesselBranch: Left(SRT, G-A101)
Mid	7.1.2.4 を参照			Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Left(SRT, G-A101)
Dist	7.1.2.4 を参照			Topographical Modifier: Distal(SRT, G-A119) VesselBranch: Left(SRT, G-A101)
Lt.Vert A	SRT	T-45700	Vertebral Artery	
	7.1.2.4 を参照			VesselBranch: Left(SRT, G-A101)

7.1.2.6 Artery of Neck (Right Extracranial Arteries, Carotid Ratios)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Artery of Neck (Right Extracranial Arteries, Carotid Ratios)	SRT	T-45005	Artery of neck	Laterality: Right(SRT, G-A100)
Rt.CCA	SRT	T-45100	Common Carotid Artery	
Prox	7.1.2.4 を参照			Topographical Modifier: Proximal(SRT, G-A118) VesselBranch: Right(SRT, G-A100)
Prox ICA/CCA(S)	LN	33868-1	ICA/CCA velocity ratio	Topographical Modifier: Proximal(SRT, G-A118)
Prox ICA/CCA(D)	99KMHC	KHU-103-00001	ICA/CCA velocity ratio(Ved)	Topographical Modifier: Proximal(SRT, G-A118)
Mid	7.1.2.4 を参照			Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Right(SRT, G-A100)
Mid ICA/CCA(S)	LN	33868-1	ICA/CCA velocity ratio	Topographical Modifier: Mid-longitudinal(SRT, G-A188)
Mid ICA/CCA(D)	99KMHC	KHU-103-00001	ICA/CCA velocity ratio(Ved)	Topographical Modifier: Mid-longitudinal(SRT, G-A188)
Dist	7.1.2.4 を参照			Topographical Modifier: Distal(SRT, G-A119) VesselBranch: Right(SRT, G-A100)
Dist ICA/CCA(S)	LN	33868-1	ICA/CCA velocity ratio	Topographical Modifier: Distal(SRT, G-A119)
Dist ICA/CCA(D)	99KMHC	KHU-103-00001	ICA/CCA velocity ratio(Ved)	Topographical Modifier: Distal(SRT, G-A119)

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7.1.2.7 Artery of Abdomen (Unilateral Abdominal Arteries (Unilateral))

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Artery of Abdomen (Unilateral Abdominal Arteries (Unilateral))	SRT	T-46002	Artery of Abdomen	Laterality: Unilateral(SRT, G-A103)
Aorta Prox	SRT	T-42000	Aorta	
Prox Ao O	SRT	G-0365	Vessel Outside Diameter	Topographical Modifier: Proximal(SRT, G-A118) VesselBranch: Unilateral(SRT, G-A103)
Prox Ao L	SRT	G-0364	Vessel Lumen Diameter	Topographical Modifier: Proximal(SRT, G-A118) VesselBranch: Unilateral(SRT, G-A103)
Prox Ao I	SRT	R-1025C	Vessel Intimal Diameter	Topographical Modifier: Proximal(SRT, G-A118) VesselBranch: Unilateral(SRT, G-A103)
Aorta Distal	SRT	T-42000	Aorta	
Distal Ao O	SRT	G-0365	Vessel Outside Diameter	Topographical Modifier: Distal(SRT, G-A119) VesselBranch: Unilateral(SRT, G-A103)
Distal Ao L	SRT	G-0364	Vessel Lumen Diameter	Topographical Modifier: Distal(SRT, G-A119) VesselBranch: Unilateral(SRT, G-A103)
Distal Ao I	SRT	R-1025C	Vessel Intimal Diameter	Topographical Modifier: Distal(SRT, G-A119) VesselBranch: Unilateral(SRT, G-A103)

Aorta Mid Suprarenal	SRT	T-42510	Supra-renal Aorta	
Mid-S Ao O	SRT	G-0365	Vessel Outside Diameter	Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Unilateral(SRT, G-A103)
Mid-S Ao L	SRT	G-0364	Vessel Lumen Diameter	Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Unilateral(SRT, G-A103)
Mid-S Ao I	SRT	R-1025C	Vessel Intimal Diameter	Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Unilateral(SRT, G-A103)
Aorta Mid Infrarenal	SRT	T-42520	Infra-renal Aorta	
Mid-I Ao O	SRT	G-0365	Vessel Outside Diameter	Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Unilateral(SRT, G-A103)
Mid-I Ao L	SRT	G-0364	Vessel Lumen Diameter	Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Unilateral(SRT, G-A103)
Mid-I Ao I	SRT	R-1025C	Vessel Intimal Diameter	Topographical Modifier: Mid-longitudinal(SRT, G-A188) VesselBranch: Unilateral(SRT, G-A103)
Aorta BIF	99KMHC	KHU-500-00001	Aortic Bifurcation	
Bif. Ao O	SRT	G-0365	Vessel Outside Diameter	VesselBranch: Unilateral(SRT, G-A103)
Bif. Ao L	SRT	G-0364	Vessel Lumen Diameter	VesselBranch: Unilateral(SRT, G-A103)
Bif. Ao I	SRT	R-1025C	Vessel Intimal Diameter	VesselBranch: Unilateral(SRT, G-A103)

7.1.2.8 Artery of Lower Extremity (Left Lower Extremity Arteries)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Artery of Lower Extremity (Left Lower Extremity Arteries)	SRT	T-47040	Artery of Lower Extremity	Laterality: Left(SRT, G-A101)
Iliac Diam	SRT	T-46710	Common Iliac Artery	
Lt.Iliac O	SRT	G-0365	Vessel Outside Diameter	VesselBranch: Left(SRT, G-A101)
Lt.Iliac L	SRT	G-0364	Vessel Lumen Diameter	VesselBranch: Left(SRT, G-A101)
Lt.Iliac I	SRT	R-1025C	Vessel Intimal Diameter	VesselBranch: Left(SRT, G-A101)

7.1.2.9 Artery of Lower Extremity (Right Lower Extremity Arteries)

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Artery of Lower Extremity (Right Lower Extremity Arteries)	SRT	T-47040	Artery of Lower Extremity	Laterality: Right(SRT, G-A100)
Iliac Diam	SRT	T-46710	Common Iliac Artery	
Rt.Iliac O	SRT	G-0365	Vessel Outside Diameter	VesselBranch: Right(SRT, G-A100)
Rt.Iliac L	SRT	G-0364	Vessel Lumen Diameter	VesselBranch: Right(SRT, G-A100)
Rt.Iliac I	SRT	R-1025C	Vessel Intimal Diameter	VesselBranch: Right(SRT, G-A100)

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7.1.3.3 Echocardiography Patient Characteristics

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Age	DCM	121033	Subject Age	
Gender	DCM	121032	Subject Sex	
BSA	LN	8277-6	Body Surface Area	
BSA Equation,DuBois	DCM	122241	$BSA = 0.007184 \cdot WT^0.725$ $0.425 \cdot HT^0.725$	

7.1.3.4 Left Ventricle

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
EDVbp	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
ESVbp	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
SVbp	SRT	F-32120	Stroke Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
CObp	SRT	F-32100	Cardiac Output	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
EFbp	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
SIbp	SRT	F-00078	Stroke Index	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
CIbp	SRT	F-32110	Cardiac Index	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
HR	LN	8867-4	Heart Rate	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2)
LVLdDiff	99KMHC	KHU-100-00001	Difference of Left Ventricular Length at End-Diastole	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVAAd2	SRT	G-0375	Left Ventricular Diastolic Area	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVAAs2	SRT	G-0374	Left Ventricular Systolic Area	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B) Cardiac Cycle Point:End Systole(SRT,R-FAB5B)
EDVAp2	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B)
ESVAp2	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B)
SVAp2	SRT	F-32120	Stroke Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B)
COAp2	SRT	F-32100	Cardiac Output	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B)

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EFap2	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B)
Slap2	SRT	F-00078	Stroke Index	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B)
Clap2	SRT	F-32110	Cardiac Index	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B)
LVA4	SRT	G-0375	Left Ventricular Diastolic Area	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVA4s	SRT	G-0374	Left Ventricular Systolic Area	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:End Systole(SRT,R-FAB5B)
EDVap4	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
ESVap4	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
SVap4	SRT	F-32120	Stroke Volume	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
COap4	SRT	F-32100	Cardiac Output	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
EFap4	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
Slap4	SRT	F-00078	Stroke Index	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
Clap4	SRT	F-32110	Cardiac Index	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
LVLd2	LN	18077-8	Left Ventricle diastolic major axis	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVLs2	LN	18076-0	Left Ventricle systolic major axis	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Two Chamber(SRT,G-A19B) Cardiac Cycle Point:End Systole(SRT,R-FAB5B)

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LVLd4	LN	18077-8	Left Ventricle diastolic major axis	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVLs4	LN	18076-0	Left Ventricle systolic major axis	Measurement Method:Method of Disks, Biplane(DCM,125207) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:End Systole(SRT,R-FAB5B)
LVLd	LN	18077-8	Left Ventricle diastolic major axis	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVAAd	SRT	G-0375	Left Ventricular Diastolic Area	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVLs	LN	18076-0	Left Ventricle systolic major axis	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:Peak Systolic(SRT,F-32021)
LVAAs	SRT	G-0374	Left Ventricular Systolic Area	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C) Cardiac Cycle Point:Peak Systolic(SRT,F-32021)
EDVap	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
ESVap	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
SVap	SRT	F-32120	Stroke Volume	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
COap	SRT	F-32100	Cardiac Output	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
EFap	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
Slap	SRT	F-00078	Stroke Index	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
Clap	SRT	F-32110	Cardiac Index	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2) Image View:Apical Four Chamber(SRT,G-A19C)
HR	LN	8867-4	Heart Rate	Measurement Method:Method of Disks, Single Plane(DCM,125208) Image Mode:2D mode(SRT,G-03A2)

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LVLd apical	LN	18077-8	Left Ventricle diastolic major axis	Image Mode:2D mode(SRT,G-03A2) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LvdsaxMv	SRT	G-0375	Left Ventricular Diastolic Area	Image Mode:2D mode(SRT,G-03A2) Image View:Parasternal short axis at the Mitral Valve level(SRT,G-039A) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LvdsaxPm	SRT	G-0375	Left Ventricular Diastolic Area	Image Mode:2D mode(SRT,G-03A2) Image View:Parasternal short axis at the Papillary Muscle level(SRT,G-039B) Cardiac Cycle Point:End Diastole(SRT,F-32011)
LVLs apical	LN	18076-0	Left Ventricle systolic major axis	Image Mode:2D mode(SRT,G-03A2) Cardiac Cycle Point:Peak Systolic(SRT,F-32021)
LvssaxMv	SRT	G-0374	Left Ventricular Systolic Area	Image Mode:2D mode(SRT,G-03A2) Image View:Parasternal short axis at the Mitral Valve level(SRT,G-039A) Cardiac Cycle Point:Peak Systolic(SRT,F-32021)
LvssaxPm	SRT	G-0374	Left Ventricular Systolic Area	Image Mode:2D mode(SRT,G-03A2) Image View:Parasternal short axis at the Papillary Muscle level(SRT,G-039B) Cardiac Cycle Point:Peak Systolic(SRT,F-32021)
EDV	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
ESV	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
SV	SRT	F-32120	Stroke Volume	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
CO	SRT	F-32100	Cardiac Output	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
EF	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
SI	SRT	F-00078	Stroke Index	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
CI	SRT	F-32110	Cardiac Index	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
HR	LN	8867-4	Heart Rate	Measurement Method:Modified Simpson(DCM,125227) Image Mode:2D mode(SRT,G-03A2)
EDV	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)
ESV	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)
SV	SRT	F-32120	Stroke Volume	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)
CO	SRT	F-32100	Cardiac Output	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)

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EF	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)
SI	SRT	F-00078	Stroke Index	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)
CI	SRT	F-32110	Cardiac Index	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)
HR	LN	8867-4	Heart Rate	Measurement Method:Bullet Method(DCM,125228) Image Mode:2D mode(SRT,G-03A2)
IVSd:2D	LN	18154-5	Interventricular Septum Diastolic Thickness	Image Mode:2D mode(SRT,G-03A2)
LVIDd:2D	LN	29436-3	Left Ventricle Internal End Diastolic Dimension	Image Mode:2D mode(SRT,G-03A2)
LVPWd:2D	LN	18152-9	Left Ventricle Posterior Wall Diastolic Thickness	Image Mode:2D mode(SRT,G-03A2)
IVSs:2D	LN	18158-6	Interventricular Septum Systolic Thickness	Image Mode:2D mode(SRT,G-03A2)
LVIDs:2D	LN	29438-9	Left Ventricle Internal Systolic Dimension	Image Mode:2D mode(SRT,G-03A2)
LVPWs:2D	LN	18156-0	Left Ventricle Posterior Wall Systolic Thickness	Image Mode:2D mode(SRT,G-03A2)
EDV:2D	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
ESV:2D	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
SV:2D	SRT	F-32120	Stroke Volume	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
CO:2D	SRT	F-32100	Cardiac Output	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
EF:2D	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
FS:2D	LN	18051-3	Left Ventricular Fractional Shortening	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
SI:2D	SRT	F-00078	Stroke Index	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
CI:2D	SRT	F-32110	Cardiac Index	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
HR:2D	LN	8867-4	Heart Rate	Measurement Method:Cube Method(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
EDV:2D	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
ESV:2D	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)

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SV:2D	SRT	F-32120	Stroke Volume	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
CO:2D	SRT	F-32100	Cardiac Output	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
EF:2D	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
FS:2D	LN	18051-3	Left Ventricular Fractional Shortening	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
SI:2D	SRT	F-00078	Stroke Index	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
CI:2D	SRT	F-32110	Cardiac Index	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
HR:2D	LN	8867-4	Heart Rate	Measurement Method:Teichholz(DCM,125209) Image Mode:2D mode(SRT,G-03A2)
EDV:2D	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
ESV:2D	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
SV:2D	SRT	F-32120	Stroke Volume	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
CO:2D	SRT	F-32100	Cardiac Output	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
EF:2D	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
FS:2D	LN	18051-3	Left Ventricular Fractional Shortening	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
SI:2D	SRT	F-00078	Stroke Index	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
CI:2D	SRT	F-32110	Cardiac Index	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
HR:2D	LN	8867-4	Heart Rate	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:2D mode(SRT,G-03A2)
IVSd:M	LN	18154-5	Interventricular Septum Diastolic Thickness	Image Mode:M mode(SRT,G-0394)
LVIDd:M	LN	29436-3	Left Ventricle Internal End Diastolic Dimension	Image Mode:M mode(SRT,G-0394)
LVPWd:M	LN	18152-9	Left Ventricle Posterior Wall Diastolic Thickness	Image Mode:M mode(SRT,G-0394)
IVSs:M	LN	18158-6	Interventricular Septum Systolic Thickness	Image Mode:M mode(SRT,G-0394)
LVIDs:M	LN	29438-9	Left Ventricle Internal Systolic Dimension	Image Mode:M mode(SRT,G-0394)

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HR:M	LN	8867-4	Heart Rate	Measurement Method:Teichholz(DCM,125209) Image Mode:M mode(SRT,G-0394)
mVcf:M	LN	59117-2	Mean Velocity of Circumferential Fiber Shortening (Mean VcFv)	Measurement Method:Teichholz(DCM,125209) Image Mode:M mode(SRT,G-0394)
EDV:M	LN	18026-5	Left Ventricular End Diastolic Volume	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
ESV:M	LN	18148-7	Left Ventricular End Systolic Volume	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
SV:M	SRT	F-32120	Stroke Volume	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
CO:M	SRT	F-32100	Cardiac Output	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
EF:M	LN	18043-0	Left Ventricular Ejection Fraction	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
FS:M	LN	18051-3	Left Ventricular Fractional Shortening	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
SI:M	SRT	F-00078	Stroke Index	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
CI:M	SRT	F-32110	Cardiac Index	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
HR:M	LN	8867-4	Heart Rate	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
mVcf:M	LN	59117-2	Mean Velocity of Circumferential Fiber Shortening (Mean VcFv)	Measurement Method:Gibson Method(99KMHC,KHU-400-00001) Image Mode:M mode(SRT,G-0394)
IVRT	LN	59083-6	Isovolumic Relaxation Time	Image Mode:Doppler Pulsed(SRT,R-409E4)
LVET:D	LN	20222-6	Ejection Time	Image Mode:Doppler Pulsed(SRT,R-409E4)
LVOT VTI	LN	20354-7	Velocity Time Integral	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
LVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	Image Mode:2D mode(SRT,G-03A2)
LVPEP:M	LN	59085-1	Pre-Ejection Period	Image Mode:M mode(SRT,G-0394)
LVSTI:M	LN	59088-5	Pre-Ejection Period/Ejection Time Ratio	Image Mode:M mode(SRT,G-0394)
LVOT Vmax	LN	11726-7	Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
LVOT PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)

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LVOT PGmean	LN	20256-4	Mean Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
LVOT SV	SRT	F-32120	Stroke Volume	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
LVOT CO	SRT	F-32100	Cardiac Output	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
LVOT SI	SRT	F-00078	Stroke Index	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
LVOT CI	SRT	F-32110	Cardiac Index	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
LVOT HR	LN	8867-4	Heart Rate	Image Mode:Doppler Pulsed(SRT,R-409E4)
E'(m)	SRT	G-037A	Left Ventricular Peak Early Diastolic Tissue Velocity	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
A'(m)	SRT	G-037C	LV Peak Diastolic Tissue Velocity During Atrial Systole	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
E'/A'(m)	LN	59129-7	Left Ventricle E to A Tissue Velocity Ratio	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
E/E'(m)	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave	
S'(m)	SRT	G-037D	Left Ventricular Peak Systolic Tissue Velocity	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
E'(l)	SRT	G-037A	Left Ventricular Peak Early Diastolic Tissue Velocity	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
A'(l)	SRT	G-037C	LV Peak Diastolic Tissue Velocity During Atrial Systole	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
E'/A'(l)	LN	59129-7	Left Ventricle E to A Tissue Velocity Ratio	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
E/E'(l)	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave	
S'(l)	SRT	G-037D	Left Ventricular Peak Systolic Tissue Velocity	Image Mode:Tissue Doppler Imaging(SRT,P5-B0128)
Mass:2D	LN	18087-7	Left Ventricle Mass	Measurement Method:CubeMethod(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
Mass_c:2D	99KMHC	KHU-100-00011	Left Ventricular Mass (c)	Measurement Method:CubeMethod(DCM,125206) Image Mode:2D mode(SRT,G-03A2)

Mass-l:2D	99KMHC	KHU-100-00012	Left Ventricular Mass Index	Measurement Method:CubeMethod(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
Mass-l_c:2D	99KMHC	KHU-100-00013	Left Ventricular Mass Index (c)	Measurement Method:CubeMethod(DCM,125206) Image Mode:2D mode(SRT,G-03A2)
Mass:M	LN	18087-7	Left Ventricle Mass	Measurement Method:CubeMethod(DCM,125206) Image Mode:M mode(SRT,G-0394)
Mass_c:M	99KMHC	KHU-100-00011	Left Ventricular Mass (c)	Measurement Method:CubeMethod(DCM,125206) Image Mode:M mode(SRT,G-0394)
Mass-l:M	99KMHC	KHU-100-00012	Left Ventricular Mass Index	Measurement Method:CubeMethod(DCM,125206) Image Mode:M mode(SRT,G-0394)
Mass-l_c:M	99KMHC	KHU-100-00013	Left Ventricular Mass Index (c)	Measurement Method:CubeMethod(DCM,125206) Image Mode:M mode(SRT,G-0394)

7.1.3.5 Right Ventricle

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
RVAWd:2D	LN	18153-7	Right Ventricle Anterior Wall Diastolic Thickness	Image Mode:2D mode(SRT,G-03A2)
RVDd:2D	LN	20304-2	Right Ventricular Internal Diastolic Dimension	Image Mode:2D mode(SRT,G-03A2)
RVDd:M	LN	20304-2	Right Ventricular Internal Diastolic Dimension	Image Mode:M mode(SRT,G-0394)
RV diam:2D	LN	20304-2	Right Ventricular Internal Diastolic Dimension	Image Mode:2D mode(SRT,G-03A2)
RV diam:M	LN	20304-2	Right Ventricular Internal Diastolic Dimension	Image Mode:M mode(SRT,G-0394)
RVET:D	LN	20222-6	Ejection Time	Image Mode:Doppler Pulsed(SRT,R-409E4)
TR RVSP	SRT	G-0380	Right Ventricular Peak Systolic Pressure	Image Mode:Doppler Pulsed(SRT,R-409E4)

7.1.3.6 Left Atrium

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
LA diam:2D	LN	29469-4	Left Atrium Antero-posterior Systolic Dimension	Image Mode:2D mode(SRT,G-03A2)
LA diam:M	LN	29469-4	Left Atrium Antero-posterior Systolic Dimension	Image Mode:M mode(SRT,G-0394)

7.1.3.7 Right Atrium

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
TR RAP	LN	18070-3	Right Atrium Systolic Pressure	Image Mode:Doppler Pulsed(SRT,R-409E4)

7.1.3.8 Aortic Valve

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
ACS:M	LN	17996-0	Aortic Valve Cusp Separation	Image Mode:M mode(SRT,G-0394)
AO/LA:2D	99KMHC	KHU-100-00008	Aortic Root to Left Atrium Ratio	Image Mode:2D mode(SRT,G-03A2)
LA/AO:2D	LN	17985-3	Left Atrium to Aortic Root Ratio	Image Mode:2D mode(SRT,G-03A2)
AO/LA:M	99KMHC	KHU-100-00008	Aortic Root to Left Atrium Ratio	Image Mode:M mode(SRT,G-0394)
LA/AO:M	LN	17985-3	Left Atrium to Aortic Root Ratio	Image Mode:M mode(SRT,G-0394)

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AV VTI	LN	20354-7	Velocity Time Integral	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
AV Vmax	LN	11726-7	Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
AV Vmean	LN	20352-1	Mean Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
AV PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
AV PGmean	LN	20256-4	Mean Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
AVA(VTI)	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Continuity Equation by Velocity Time Integral(DCM,125215) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
AVA(Vmax)	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Continuity Equation by Peak Velocity(DCM,125214) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
AVA(Trace)	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Planimetry(DCM,125220) Image Mode:Doppler Pulsed(SRT,R-409E4)
LVPEP:D	LN	59085-1	Pre-Ejection Period	Image Mode:Doppler Pulsed(SRT,R-409E4)
LVAccTime	LN	20168-1	Acceleration Time	Image Mode:Doppler Pulsed(SRT,R-409E4)
LV PEP/ET	LN	59088-5	Pre-Ejection Period/Ejection Time Ratio	Image Mode:Doppler Pulsed(SRT,R-409E4)
LV Acce T/ET	SRT	G-0382	Ratio of Aortic Valve Acceleration Time to Ejection Time	Image Mode:Doppler Pulsed(SRT,R-409E4)
Decel Rate	LN	20216-8	Deceleration Slope	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
AR Vmax	LN	11726-7	Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
AR PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
AR V ed	LN	11653-3	End Diastolic Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
AR Time	LN	20217-6	Deceleration Time	Image Mode:Doppler Pulsed(SRT,R-409E4)
AI Vmax	99KMHC	KHU-100-00014	Aortic Insufficiency maximum Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
AI PGmax	99KMHC	KHU-100-00015	Aortic Insufficiency maximum Pressure Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)

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AI PHT	LN	20280-4	Pressure Half-Time	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
AI Dec Time	99KMHC	KHU-100-00016	Aortic Insufficiency Deceleration Time	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
AI Dec Slope	99KMHC	KHU-100-00017	Aortic Insufficiency Deceleration Slope	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)

7.1.3.9 Mitral Valve

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
EPSS:2D	LN	59098-4	Mitral Valve E-septal Separation	Image Mode:2D mode(SRT,G-03A2)
MV E pt	LN	18037-2	Mitral Valve E-Wave Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4)
MV A pt	LN	17978-8	Mitral Valve A-Wave Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4)
Dec Time	SRT	G-0384	Mitral Valve E-Wave Deceleration Time	Image Mode:Doppler Pulsed(SRT,R-409E4)
Dec Slope	LN	20216-8	Deceleration Slope	Image Mode:Doppler Pulsed(SRT,R-409E4)
MV Vmax	LN	11726-7	Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MV Vmean	LN	20352-1	Mean Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MV PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MV PGmean	LN	20256-4	Mean Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MV VTI	LN	20354-7	Velocity Time Integral	Image Mode:Doppler Pulsed(SRT,R-409E4)
MV PHT	LN	20280-4	Pressure Half-Time	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MV diam	SRT	G-038F	Cardiovascular Orifice Diameter	Image Mode:Doppler Pulsed(SRT,R-409E4)
MV E duration	99KMHC	KHU-100-00002	Mitral Valve E-Wave Duration	Image Mode:Doppler Pulsed(SRT,R-409E4)
MV A duration	SRT	G-0385	Mitral Valve A-Wave Duration	Image Mode:Doppler Pulsed(SRT,R-409E4)

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MV SV	SRT	F-32120	Stroke Volume	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
MV CO	SRT	F-32100	Cardiac Output	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
MV SI	SRT	F-00078	Stroke Index	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
MV CI	SRT	F-32110	Cardiac Index	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)
MV E/A:D	LN	18038-0	Mitral Valve E to A Ratio	Image Mode:Doppler Pulsed(SRT,R-409E4)
MV A/E:D	99KMHC	KHU-100-00003	Mitral Valve A to E Ratio	Image Mode:Doppler Pulsed(SRT,R-409E4)
MVA(PHT)	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Area by PHT(DCM,125210) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MVA(VTI)	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Continuity Equation(DCM,125212) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MVA(Trace)	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Planimetry(DCM,125220) Image Mode:2D mode(SRT,G-03A2) Direction of Flow:Antegrade Flow(SRT,R-42047)
MV HR	LN	8867-4	Heart Rate	Image Mode:Doppler Pulsed(SRT,R-409E4)
E/E'(m)	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave	
E/E'(l)	SRT	G-037B	Ratio of MV Peak Velocity to LV Peak Tissue Velocity E-Wave	
MV E amp:M	99KMHC	KHU-100-00004	Mitral Valve E-Wave Amplitude	Image Mode:M mode(SRT,G-0394)
MV A amp:M	99KMHC	KHU-100-00005	Mitral Valve A-Wave Amplitude	Image Mode:M mode(SRT,G-0394)
DEexcursion	LN	17997-8	Mitral Valve D-E Excursion	Image Mode:M mode(SRT,G-0394)
DE amp:M	99KMHC	KHU-100-00006	D-E Amplitude of the Anterior Mitral Leaflet	Image Mode:M mode(SRT,G-0394)
EPSS:M	LN	18036-4	Mitral Valve EPSS, E wave	Image Mode:M mode(SRT,G-0394)
EF slope:M	LN	18040-6	Mitral Valve E-F Slope by M-Mode	Image Mode:M mode(SRT,G-0394)
MV A/E amp:M	99KMHC	KHU-100-00007	Mitral Valve ratio of the A point to the E point	Image Mode:M mode(SRT,G-0394)

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MR Vmax	LN	11726-7	Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
MR PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
dt	LN	20217-6	Deceleration Time	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
dP/dt	LN	18035-6	Mitral Regurgitation dP/dt derived from Mitral Reg. velocity	Image Mode:Doppler Pulsed(SRT,R-409E4)
MR Radius	LN	59102-4	Flow Radius	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
MRAliasVel	LN	59130-5	Alias Velocity	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
MR VTI	LN	20354-7	Velocity Time Integral	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
MRFlowRate	LN	34141-2	Peak Instantaneous Flow Rate	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
EO Area	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
MR Flow Vol	LN	33878-0	Volume Flow	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
MS Radius	LN	59102-4	Flow Radius	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MSAliasVel	LN	59130-5	Alias Velocity	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MS VTI	LN	20354-7	Velocity Time Integral	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)

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MS Vmax	LN	11726-7	Peak Velocity	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MS PGmax	LN	20247-3	Peak Gradient	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MS Angle	99KMHC	KHU-100-00010	Angle	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MSFlowRate	LN	34141-2	Peak Instantaneous Flow Rate	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MS Flow Vol	LN	33878-0	Volume Flow	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
MVA	SRT	G-038E	Cardiovascular Orifice Area	Measurement Method:Proximal Isovelocity Surface Area(DCM,125216) Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)

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RVOT Vmean	LN	20352-1	Mean Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT PGmean	LN	20256-4	Mean Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT diam	SRT	G-038F	Cardiovascular Orifice Diameter	Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT SV	SRT	F-32120	Stroke Volume	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT CO	SRT	F-32100	Cardiac Output	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT SI	SRT	F-00078	Stroke Index	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT CI	SRT	F-32110	Cardiac Index	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4) Image View:Right Ventricular Outflow Tract View(SRT,G-039D)
RVOT HR	LN	8867-4	Heart Rate	Image Mode:Doppler Pulsed(SRT,R-409E4)
RVPEP:D	LN	20301-8	Right Ventricle Pre-Ejection Period	Image Mode:Doppler Pulsed(SRT,R-409E4)
RVAccTime	LN	20168-1	Acceleration Time	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
RV PEP/ET	LN	59088-5	Pre-Ejection Period/Ejection Time Ratio	Image Mode:Doppler Pulsed(SRT,R-409E4)
RV Acce T/ET	SRT	G-0388	Ratio of Pulmonic Valve Acceleration Time to Ejection Time	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,R-42E61)
PREDV	LN	11653-3	End Diastolic Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:End Diastole(SRT,F-32011)
PREDPG	99KMHC	KHU-100-00009	End Diastolic Pressure Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:End Diastole(SRT,F-32011)

7.1.3.11 Tricuspid Valve

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
TV VTI	LN	20354-7	Velocity Time Integral	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
TV Vmax	LN	11726-7	Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
TV PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
TV Vmean	LN	20352-1	Mean Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
TV PGmean	LN	20256-4	Mean Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Antegrade Flow(SRT,R-42047)
TV E pt	LN	18031-5	Tricuspid Valve E Wave Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4)
TV A pt	LN	18030-7	Tricuspid Valve A Wave Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4)
TV E/A	LN	18039-8	Tricuspid Valve E to A Ratio	Image Mode:Doppler Pulsed(SRT,R-409E4)
TR Vmax	LN	11726-7	Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
TR Vmean	LN	20352-1	Mean Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
TR PGmax	LN	20247-3	Peak Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
TR PGmean	LN	20256-4	Mean Gradient	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)
TR VTI	LN	20354-7	Velocity Time Integral	Image Mode:Doppler Pulsed(SRT,R-409E4) Direction of Flow:Regurgitant Flow(SRT,G-0367)

7.1.3.12 Aorta

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
AO:2D	LN	18015-8	Aortic Root Diameter	Image Mode:2D mode(SRT,G-03A2)
AO:M	LN	18015-8	Aortic Root Diameter	Image Mode:M mode(SRT,G-0394)

7.1.3.13 Pulmonary Venous Structure

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
PVs1 Vel	99KMHC	KHU-100-00018	Pulmonary Vein First Systolic Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Systolic(SRT,F-32020)
PVs2 Vel	LN	29450-4	Pulmonary Vein Systolic Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Systolic(SRT,F-32020)
PVd Vel	LN	29451-2	Pulmonary Vein Diastolic Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Diastole(SRT,F-32010)
PVa Vel	LN	29453-8	Pulmonary Vein Atrial Contraction Reversal Peak Velocity	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Diastole(SRT,F-32010)
PVa duration	SRT	G-038B	Pulmonary Vein A-Wave Duration	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Atrial Systole(SRT,F-32030)
PVs VTI	SRT	G-038C	Pulmonary Vein S-Wave Velocity Time Integral	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Systolic(SRT,F-32020)
PVd VTI	SRT	G-038D	Pulmonary Vein D-Wave Velocity Time Integral	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Diastole(SRT,F-32010)
PVs2/PVd	LN	29452-0	Pulmonary Vein Systolic to Diastolic Ratio	Image Mode:Doppler Pulsed(SRT,R-409E4)
PVd DecT	LN	20217-6	Deceleration Time	Image Mode:Doppler Pulsed(SRT,R-409E4) Cardiac Cycle Point:Diastole(SRT,F-32010)
Sys Fraction	LN	59113-1	Pulmonary Vein A VTI to Mitral Valve VTI Ratio	Image Mode:Doppler Pulsed(SRT,R-409E4)

7.1.3.14 Cardiac Shunt Study

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
Qp/Qs	LN	29462-9	Pulmonary-to-Systemic Shunt Flow Ratio	Measurement Method:Doppler Volume Flow(DCM,125219) Image Mode:Doppler Pulsed(SRT,R-409E4)

7.1.3.15 Vena Cava

Label	Coding Scheme Designator (0008,0102)	Code Value (0008,0100)	Code Meaning (0008,0104)	Modifiers
E-diam:2D	LN	18006-7	Inferior Vena Cava Diameter	Image Mode:2D mode(SRT,G-03A2) Respiratory Cycle Point:Expiration(SRT,F-20020)
I-diam:2D	LN	18006-7	Inferior Vena Cava Diameter	Image Mode:2D mode(SRT,G-03A2) Respiratory Cycle Point:Inspiration(SRT,F-20010)
E-diam:M	LN	18006-7	Inferior Vena Cava Diameter	Image Mode:M mode(SRT,G-0394) Respiratory Cycle Point:Expiration(SRT,F-20020)
I-diam:M	LN	18006-7	Inferior Vena Cava Diameter	Image Mode:M mode(SRT,G-0394) Respiratory Cycle Point:Inspiration(SRT,F-20010)

KONICA MINOLTA

A747EA01JA11

2020-09-11
(AI)